# Posey® Secure Tie For Tracheostomy Tube – 8196S, 8196M, 8196L


**DESCRIPTION OF PRODUCT:** Posey Secure Tie for Tracheostomy Tube: Three sizes fit tiny preemies to large adults. Ensure proper fit before use.

**INTENDED USE:** The Posey Secure Tie is used to hold a tracheal tube in place.

**CONTRAINDICATIONS:** None

**APPLICATION INSTRUCTIONS:**

1. Thread the “hook” tape through one slot of the tracheostomy tube and attach the “hook” fastener to the fuzzy side of the neck pad or strap (See Fig. 1). Slot should be at least 3/8” (1 cm) to allow for easy insertion.
2. Bring the foam neck collar behind the patient’s neck and thread “hook” tape through other slot in the tracheostomy tube. Secure “hook” fastener to fuzzy side of neck strap or neck collar.
3. Crisscross the white Secure Ties underneath the patient’s connector (swivel/adaptor) or flex hose that attaches to the trach tube (See Fig. 2).
4. Grasping each end of the Secure Tie, bring over the top of the patient connector (swivel adaptor/elbow) or flex hose. The Secure Tie may be tied at this point in a single tie or a bow tie around the patient connector as illustrated (See Fig. 3).

**POSEY SECURE TIE FOR TRACHEOSTOMY TUBE (One dozen per box):**

**REF 8196S** Small 7-9”, (18-23 cm)

**REF 8196M** Medium 9-14”, (23-35 cm)

**REF 8196L** Large 13-19”, (33-48 cm)

**WARNING**

- **MONITORING:** Be sure to follow your facility’s policies and guidelines for frequency of patient monitoring, for proper product application, and skin integrity. ALWAYS monitor patients carefully to ensure they are not at risk for disrupting lifesaving therapy. Constant supervision, electronic monitoring and/or the use of a restraint may be required in conjunction with proper medical supervision per your facility’s policies.
- **INSPECTION:** Inspect before each use; check for broken stitches or parts, or torn, cut or frayed material. DO NOT use soiled or damaged product.
- **REMOVAL:** Do not force tabs through the slot when removing trach tie. This could cause accidental extubation or displacement of the tube. If the slot is too small to allow for easy removal, cut off the hook and loop tab, and pull the tie gently.
- Discard if “hook” does not firmly attach to fuzzy backing.
- Dispose of the product per facility’s policy if worn or damaged.

EN Report any serious incident related to device use to TIDI Products and the Member State Competent Authority.

Translations of this document in French, Spanish, Dutch, German, Italian, and Portuguese are located at [www.tidiproducts.com/IFU](http://www.tidiproducts.com/IFU).

TIDI Products, LLC • 570 Enterprise Drive, Neenah, WI 54956 USA  
Phone: 1.800.447.6739 • International: +1.920.751.4036 • [www.tidiproducts.com](http://www.tidiproducts.com)

UK CA UK Responsible Person  
Emergo Consulting (UK) Limited  
c/o Cr360 – UL International  
Compass House, Vision Park Histon  
Cambridge CB249BZ  
United Kingdom

M1395 REV D 2022-01-26

© Posey Products, LLC. All rights reserved.

# Posey® Secure Tie For Tracheostomy Tube – 8196S, 8196M, 8196L


**DESCRIPTION OF PRODUCT:** Posey Secure Tie for Tracheostomy Tube: Three sizes fit tiny preemies to large adults. Ensure proper fit before use.

**INTENDED USE:** The Posey Secure Tie is used to hold a tracheal tube in place.

**CONTRAINDICATIONS:** None

**APPLICATION INSTRUCTIONS:**

1. Thread the “hook” tape through one slot of the tracheostomy tube and attach the “hook” fastener to the fuzzy side of the neck pad or strap (See Fig. 1). Slot should be at least 3/8” (1 cm) to allow for easy insertion.
2. Bring the foam neck collar behind the patient’s neck and thread “hook” tape through other slot in the tracheostomy tube. Secure “hook” fastener to fuzzy side of neck strap or neck collar.
3. Crisscross the white Secure Ties underneath the patient’s connector (swivel/adaptor) or flex hose that attaches to the trach tube (See Fig. 2).
4. Grasping each end of the Secure Tie, bring over the top of the patient connector (swivel adaptor/elbow) or flex hose. The Secure Tie may be tied at this point in a single tie or a bow tie around the patient connector as illustrated (See Fig. 3).

**POSEY SECURE TIE FOR TRACHEOSTOMY TUBE (One dozen per box):**

**REF 8196S** Small 7-9”, (18-23 cm)

**REF 8196M** Medium 9-14”, (23-35 cm)

**REF 8196L** Large 13-19”, (33-48 cm)

**WARNING**

- **MONITORING:** Be sure to follow your facility’s policies and guidelines for frequency of patient oring, for proper product application, and skin integrity. ALWAYS monitor patients carefully to ensure they are not at risk for disrupting lifesaving therapy. Constant supervision, electronic monitoring and/or the use of a restraint may be required in conjunction with proper medical supervision per your facility’s policies.
- **INSPECTION:** Inspect before each use; check for broken stitches or parts, or torn, cut or frayed material. DO NOT use soiled or damaged product.
- **REMOVAL:** Do not force tabs through the slot when removing trach tie. This could cause accidental extubation or displacement of the tube. If the slot is too small to allow for easy removal, cut off the hook and loop tab, and pull the tie gently.
- Discard if “hook” does not firmly attach to fuzzy backing.
- Dispose of the product per facility’s policy if worn or damaged.

EN Report any serious incident related to device use to TIDI Products and the Member State Competent Authority.

Translations of this document in French, Spanish, Dutch, German, Italian, and Portuguese are located at [www.tidiproducts.com/IFU](http://www.tidiproducts.com/IFU).

TIDI Products, LLC • 570 Enterprise Drive, Neenah, WI 54956 USA  
Phone: 1.800.447.6739 • International: +1.920.751.4036 • [www.tidiproducts.com](http://www.tidiproducts.com)

UK CA UK Responsible Person  
Emergo Consulting (UK) Limited  
c/o Cr360 – UL International  
Compass House, Vision Park Histon  
Cambridge CB249BZ  
United Kingdom

M1395 REV D 2022-01-26

© Posey Products, LLC. All rights reserved.

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TIDI Products, LLC • 570 Enterprise Drive, Neenah, WI 54956 EUA
**CE REP** MDSS GmbH
Emurgo Consulting (UK) Limited
© Posey Products, LLC. Todos os direitos reservados.