REF MCGLPICC

- (EN) Catheter Fixation for MedComp® PICC/CVC
- (DE) Katheter Fixierung für MedComp® PICC/CVC
- (FR) Fixation de Cathéter pour MedComp® PICC/CVC
- (IT) Fissaggio del catetere MedComp® PICC/CVC
- (NL) Catheter Fixatie voor MedComp® PICC/CVC
- (SV) Kateter Fixering för MedComp® PICC/CVC
- (ES) La Fijación del Catéter para el MedComp® PICC/CVC
- (PT) Fixação de catéter para MedComp® PICC/CVC
- (EL) Ηθέσμεση για MedComp® PICC/CVC
- (FI) Katetri Pakkomielle MedComp® PICC/CVC
- (DA) Kateter Fiksering for MedComp® PICC/CVC
- (S) Catheter Fixation pentru MedComp® PICC/CVC
- (NO) Kateter Fiksering for MedComp® PICC/CVC
- (PL) Catheter Utrwalania dla MedComp® PICC/CVC

Grip-Lok®

44U00012 Rev 01

2022-10-13

## INSTRUCTIONS

EN  
English

Intended for use in the general population.

1. **Select** the area for the placement of the Grip-Lok. Note: The catheter should be already inserted in the patient.
2. **Prepare** the skin according to the standard hospital protocol for dressing application. Hair removal may be required on some patients for better adhesion.
3. **Open** the top flap, remove the interior liner and hold the flap open.
4. **Slide** the Grip-Lok under the catheter and adhere in the exposed adhesive cut-out space provided.
5. **Secure** the top flap over the catheter.
6. **Hold** the Grip-Lok in position while removing the paper backing from one side, then the other, to secure on the skin surface.

**Direction:** For additional skin protection, use skin prep pad prior to application of stabilization device.**Note:** Use of an alcohol swab may help removal of securement device from the skin.

- ⚠ If you experience severe redness, itching, swelling, or irritation of the skin, consult your physician as this may be a sign of an allergic reaction.
- ⚠ Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.
- ⚠ Replace securement device if soiled or saturated in fluid or if device shows signs of wear or damage.
- Note:** If being used in a medical facility, dispose according to medical facility and municipal guidelines.
- ⊗ Do not use if the product's sterile barrier system or its packaging is compromised. Report any serious incident related to device use to TIDI Products and the Member State competent authority.

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TIDI® Products, LLC  
570 Enterprise Drive  
Neenah, WI 54956 USA  
Phone: +1 920 751 4300  
tidiproducts.com/securement

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UK Responsible Person:  
Emergo Consulting (UK) Limited  
c/o Cr 360 - UL International  
Compass House, Vision Park Histon  
Cambridge CB24 9BZ  
United Kingdom

Medical Components, Inc.  
1499 Delp Drive  
Harleysville, PA 19438  
Phone: 215.256.4201  
Fax: 215.256.1787  
www.medcompnet.com

Grip-Lok is a registered trademark of TIDI Products, LLC.

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Grip-Lok jest zastrzeżonym znakiem towarowym firmy TIDI Products, LLC.
MDSS GmbH