## INSTRUCTIONS

EN

Intended for use in the general population, for adhesive securement of the catheter hub, tube or line to the patient's skin. English

- Select** the area for the placement of the Grip-Lok. Note: The catheter should be already inserted in the patient.
- Prepare** the skin according to the standard hospital protocol for dressing application. Skin prep or hair removal may be required on some patients for better adhesion.
- Open** the top flap and hold the flap open.
- Position** the Grip-Lok with the catheter in the groove at edge of the foam section.
- Hold** the Grip-Lok in position while removing the paper backing from one side, then the other, to secure on the skin surface.
- Remove** the interior liner and place the catheter over the foam and onto the center of the exposed adhesive surface.
- Secure** the top flap over the catheter and apply gentle pressure to the center fabric area over the catheter adhesive.

**Direction:** for added skin protection, skin prep may be used prior to application of securement device.

**Note:** Use of an alcohol swab may help removal of securement device from the skin.

⚠ If you experience severe redness, itching, swelling, or irritation of the skin, consult your physician as this may be a sign of an allergic reaction.

⚠ Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

⚠ Replace securement device if soiled or saturated in fluid or if device shows signs of wear or damage.

**Note:** If being used in a medical facility, dispose according to medical facility and municipal guidelines.

⊘ Do not use if the product's sterile barrier system or its packaging is compromised.

Report any serious incident related to device use to TIDI Products and the Member State competent authority.

HYPOALLERGENIC

TIDI Products, LLC  
570 Enterprise Drive  
Neenah, WI 54956, USA  
Phone: + 1.920.751.4300  
tidiproducts.com/securement

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UK Responsible Person:  
Emergo Consulting (UK) Limited  
c/o Cr 360 - UL International  
Compass House, Vision Park Histon  
Cambridge CB24 9BZ  
United Kingdom




Grip-Lok is a registered trademark of TIDI Products, LLC.  
May be patented: see [www.tidiproducts.com/virtual-patent-marking](http://www.tidiproducts.com/virtual-patent-marking)

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