## INSTRUCTIONS

EN

Intended for use in the general population, for adhesive securement of the catheter hub, tube or line to the patient's skin. English

1. **Select** the area for the placement of the Grip-Lok. Note: The catheter should be already inserted in the patient with Luer connection attached.
2. **Prepare** the skin according to the standard hospital protocol for dressing application.
3. **Remove** one side of the release liner from the large portion of the Grip-Lok.
4. **Place** the Grip-Lok centered on the nose with the thin tab portions extending over the end of the nose.
5. **Press** the exposed adhesive down one side of the nose and onto the face.
6. **Remove** other side of the release liner from the large portion of the Grip-Lok and press the exposed adhesive down the other side of the nose and onto the face.
7. **Remove** the adhesive release liner from under one side of the thin tab portions extending over the end of the nose and wrap around the tubing.
8. **Remove** the adhesive release liner from under the other side of the thin tab portions and wrap around the tubing in the opposite direction.

**Direction:** For added skin protection, skin prep may be used prior to application of securement device.

**Note:** Use of an alcohol swab may help removal of securement device from the skin.

**Note:** You can make removal of the Grip-Lok from the tubing easier by folding the end of the thin tab portions back against itself.

⚠ If you experience severe redness, itching, swelling, or irritation of the skin, consult your physician as this may be a sign of an allergic reaction.

⚠ Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

⚠ Replace securement device if soiled or saturated in fluid or if device shows signs of wear or damage.

**Note:** If being used in a medical facility, dispose according to medical facility and municipal guidelines.

⊘ Do not use if the product's sterile barrier system or its packaging is compromised.

Report any serious incident related to device use to TIDI Products and the Member State competent authority.

HYPOALLERGENIC

2023-08-22

TIDI Products, LLC  
570 Enterprise Drive  
Neeenah, WI 54956, USA  
Phone: + 1.920.751.4300  
tidiproducts.com/securement

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UK  
CA  
0086

UK Responsible Person:  
Emergo Consulting (UK) Limited  
c/o Cr 360 - UL International  
Compass House, Vision Park Histon  
Cambridge CB24 9BZ  
United Kingdom

Grip-Lok is a registered trademark of TIDI Products, LLC.  
May be patented: see [www.tidiproducts.com/virtual-patent-marking](http://www.tidiproducts.com/virtual-patent-marking)

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