INSTRUCTIONS

EN

Intended for use in the general population, for adhesive securing of the catheter hub, tube or line to the patient's skin.

1. Select the area for the placement of the Grip-Lok. Note: catheter should be already placed/inserted.

2. Prepare the skin according to the standard hospital protocol for dressing application. Skin must be clean and dry. Skin prep or hair removal may be required on some patients for better adhesion.

3. Open the top fabric flap and remove the interior liner to expose the adhesive areas.

4. Slide the Grip-Lok under the catheter hub centering it in the exposed adhesive area with the distal end of the hub extending just beyond the Grip-Lok.

5. Secure top flap over the catheter and apply firm pressure to the fabric area over the adhesive and around and in between lumens.

6. Pull the paper backing from the bottom of one side of the Grip-Lok, then the other, to secure in desired position on the skin.

Direction: for added skin protection, skin prep may be used prior to application of securement device.

Note: Use of an alcohol swab may help removal of securement device from the skin.

⚠️ If you experience severe redness, itching, swelling, or irritation of the skin, consult your physician as this may be a sign of an allergic reaction.

⚠️ Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

⚠️ Replace securement device if soiled or saturated in fluid or if device shows signs of wear or damage.

Note: If being used in a medical facility, dispose according to medical facility and municipal guidelines.

🚫 Do not use if the product's sterile barrier system or its packaging is compromised.

Report any serious incident related to device use to TIDI Products and the Member State competent authority.

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HYPALLERGEN

TIDI Products, LLC
570 Enterprise Drive
Neenah, WI 54956, USA
Phone: +1.920.751.4300
tidiproducts.com/securement

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UK CA 0086
UK Responsible Person:
Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

Grip-Lok is a registered trademark of TIDI Products, LLC.
May be patented: see www.tidiproducts.com/virtual-patent-marking

2023-08-22

Grip-Lok®

REF 3601CVC

(EN) Medium Universal PICC and CVC Securement Device, For lines and tubes sized 14 - 40 Fr

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Grip-Lok est une marque déposée de TIDI Products, LLC.