INSTRUCTIONS

EN

Intended for use in the general population, for adhesive securement of the catheter hub, tube or line to the patient's skin. English

1. Select the correct size Grip-Lok for your application.

REF	Size	Line/Tubing Diameter
3200S	Small	4.5 - 13.5 Fr
3300M	Medium	9 - 24 Fr
3400L	Large	16 - 40 Fr

2. Select the area for the placement of the Grip-Lok. Note: The line or tube should be already inserted in the patient.

3. Prepare the skin according to the standard hospital protocol for dressing application. Skin prep or hair removal may be required on some patients for better adhesion.

4. Open the top flap, remove the interior liner and hold the flap open.

5. Slide the Grip-Lok under so that the line or tube is in the center of the exposed adhesive.

6. Secure the top flap over the line or tube and apply gentle pressure to the center fabric area over the adhesive.

7. Hold the Grip-Lok in position while removing the paper backing from one side, then the other, to secure on the skin surface.

Direction: For additional skin protection, use skin prep pad prior to application of stabilization device.

Note: Use of an alcohol swab may help removal of securement device from the skin.

- ⚠ If you experience severe redness, itching, swelling, or irritation of the skin, consult your physician as this may be a sign of an allergic reaction.

- ⚠ Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

- ⚠ Replace securement device if soiled or saturated in fluid or if device shows signs of wear or damage.

Note: If being used in a medical facility, dispose according to medical facility and municipal guidelines.

- ⊘ Do not use if the product's sterile barrier system or its packaging is compromised.

Report any serious incident related to device use to TIDI Products and the Member State competent authority.

HYPOALLERGENIC

2023-08-10

44U0009 Rev 02

TIDI Products, LLC
570 Enterprise Drive
Neenah, WI 54956, USA
Phone: + 1.920.751.4300
tidiproducts.com/securement

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UK Responsible Person:
Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom


Grip-Lok is a registered trademark of TIDI Products, LLC

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Grip-Lok jest zastrzeżonym znakiem towarowym firmy TIDI Products, LLC.
Grip-Lok je registrovaná ochranná známka společnosti TIDI Products, LLC.