## INSTRUCTIONS

EN

Intended for use in the general population, for adhesive securement of the catheter hub, tube or line to the patient's skin. English

- Select** the area for the placement of the Grip-Lok.  
Note: The arterial cannula should be already inserted in the patient.
- Prepare** the skin according to the standard hospital protocol for dressing application. Skin prep or hair removal may be required on some patients for better adhesion.
- Open** the top flap, remove the interior liner and hold the flap open.
- Slide** the Grip-Lok under the arterial cannula and adhere in the exposed adhesive cut-out space provided.
- Hold** the Grip-Lok in position while removing the paper backing from one side, then the other, to secure on the skin surface.
- Secure** the top flap over the arterial cannula with the FloSwitch portion coming through the center cut-out.

**Direction:** For additional skin protection, use skin prep pad prior to application of stabilization device.

**Note:** Use of an alcohol swab may help removal of securement device from the skin.

- ⚠ If you experience severe redness, itching, swelling, or irritation of the skin, consult your physician as this may be a sign of an allergic reaction.
  - ⚠ Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.
  - ⚠ Replace securement device if soiled or saturated in fluid or if device shows signs of wear or damage.
- Note:** If being used in a medical facility, dispose according to medical facility and municipal guidelines.
- ⊗ Do not use if the product's sterile barrier system or its packaging is compromised.
- Report any serious incident related to device use to TIDI Products and the Member State competent authority.

HYPOALLERGENIC

TIDI® Products, LLC  
570 Enterprise Drive  
Neenah, WI 54956 USA  
Phone: +1 920 751 4300  
tidiproducts.com/securement

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UK CA 0086

UK Responsible Person:  
Emergo Consulting (UK) Limited  
c/o Cr 360 - UL International  
Compass House, Vision Park Histon  
Cambridge CB24 9BZ  
United Kingdom




BD and FloSwitch trademarks are property of Becton, Dickinson and Company.  
Grip-Lok is a registered trademark of TIDI Products, LLC.  
May be patented: see [www.tidiproducts.com/virtual-patent-marking](http://www.tidiproducts.com/virtual-patent-marking)

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Grip-Lok jest zastrzeżonym znakiem towarowym firmy TIDI Products, LLC.