# Posey® Limb Holder 2534

## Application Instructions for Wrist and Ankle

**Rx ONLY**

**DESCRIPTION OF PRODUCT:** Padded limb holder. For bed and stretcher use only.


**INTENDED USE:**

- To protect patients assessed to be at risk of disrupting life-saving treatments (e.g., chronic tube pulling) or in danger of injury to themselves or to others by limiting limb movement in a bed or stretcher.
- Follow your hospital's restraint policies and procedures that are compliant with USA CMS guidelines and state laws, or other governing agencies outside the USA.


**CONTRAINDICATIONS:** None.

**APPLICATION INSTRUCTIONS:** The order of applying/removing the cuff to the patient and the straps to the bed is to be determined by the clinician. Use method a. or b. below to attach straps to frame (repeat steps 1–3 on each side).


- 1a. Triangulation process to restrict patient's range of motion: Separate the straps and attach at different points along the frame that moves with the patient (Fig. 1), out of the patient's reach, using quick-release buckles (do not attach to side rail or head/footboard) (Fig. 2).


- 1b. To increase patient's range of motion: Place the straps together and attach to a single point along a part of the frame, out of the patient's reach, using quick-release buckle.

2. Attach the female end of the quick-release buckle (short strap) to the frame that moves with the patient (Fig. 1), out of the patient's reach (do not attach to side rail or head/footboard). You may also wrap the connecting strap once around the frame to move the buckle out of the patient's reach. Secure by feeding the female end through the loop in the strap.


3. Insert the male end of the connecting strap into the female end of the short strap. Listen for a snapping sound. Pull firmly on the straps to ensure a secure connection. To limit unwanted adjustment, tie an overhand knot with the excess strap directly below the quick-release buckle.


**APPLYING THE CUFFS:** (Repeat steps 1–4)

1. Wrap the limb holder cuff around the patient's wrist/ankle so that the buckle and connecting strap are on the ulnar side of the wrist or lateral malleolus of the ankle.

2. Secure the hook-and-loop fastener. Slide ONE finger (flat) between the cuff and the inside of the patient's wrist/ankle to ensure proper fit (Fig. 3). The strap must be snug, but not restrict circulation.

3. Close the quick-release buckle on the cuff. Insert ONE finger (flat) under the buckle and pull the strap snug, but not so tight as to restrict circulation (Fig. 4). Release the quick-release buckle, twist buckle 180°, and reconnect (Fig. 5). Listen for a snapping sound.

4. Attach the hook end of the cuff strap to the fuzzy backing on the cuff to keep the quick-release buckle from sliding (Fig. 6).

**POSEY LIMB HOLDER**

**REF 2534** Foam Quick-Release Buckle Cuff (double strap with 2-piece quick-release buckle)

**ADJUSTING STRAP LENGTH:** Adjust the connecting strap(s) length from the frame to allow desired freedom of movement without compromising patient or caregiver safety. To limit unwanted adjustment, tie an overhand knot with the excess strap directly below the quick-release buckle. Ensure the remaining strap end(s) are secured and out of the patient's reach.

**TO REMOVE:** Release the cuff from the patient. Release the anchor strap(s) from the bed frame.

**TO LIMIT LOWER LIMB RANGE OF MOTION FOR LEGS (Fig. 7):** The clinician applying the strap's criss-cross application method can determine which is applied first (cuff or strap) and which end of the restraint is criss-crossed.

1. Attach the cuff secured to the bottom right side of the frame that moves with the patient, out of the patient's reach, to the left ankle.
2. Criss-cross the straps and attach the cuff secured to the bottom left side of the frame that moves with the patient, out of the patient's reach, to the right ankle.
3. Adjust connecting straps as necessary.

**PRECAUTIONS:**

- A clinical assessment and decision are required when used with monitoring lines or if the patient has a wound or dislocated/fractured limb.
- Check the patient regularly to ensure that circulation is not impaired. Serious injury may occur if the cuff restricts circulation when the limb holder is applied.

**WARNING:** Additional or different body or limb restraints may be needed (visit [tidiproductions.com](http://tidiproductions.com))

- Some patients may require additional interventions in conjunction with a restraint in order to prevent injury to self or others:
  - If the patient pulls violently against the bed straps.
  - To reduce the risk of the patient getting access to the line/wound/tube site.
  - To prevent the patient from flailing or bucking up and down, causing self-injury.
- It is a clinical decision to decide when a restraint is no longer deemed safe.

**BED SAFETY:** Refer to the Food and Drug Administration (FDA), or other governing agencies outside the USA, for the most recent Hospital Bed Safety Guidelines and the bed manufacturer's instructions for use.

EN

Report any serious incident related to device use to TIDI Products and the Member State Competent Authority. | Translations of this document in French, Spanish, Dutch, German, Italian, and Portuguese are located at [www.tidiproductions.com/IFU](http://www.tidiproductions.com/IFU).

**ADDITIONAL SAFETY AND LAUNDERING INSTRUCTIONS ON OTHER SIDE.**

TIDI Products, LLC • 570 Enterprise Drive, Neenah, WI 54956 USA  
Phone: 1.800.521.1314 • International: +1.920.751.4036 • [www.tidiproductions.com](http://www.tidiproductions.com)

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UK Responsible Person  
Emergo Consulting (UK) Limited  
c/o CR360 – UL International  
Compass House, Vision Park, Histon,  
Cambridge CB24 9BZ  
United Kingdom

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19291 REV D 03-25-2025


**For RX PRODUCTS ONLY.** Federal law (USA) restricts RX devices to sale by or on order of a physician. For use in a licensed healthcare facility only.

**STAFF TRAINING:** TIDI offers in-service training aids at no charge. Contact TIDI online at [www.tidiproducts.com](http://www.tidiproducts.com) or call toll-free at 1.800.521.1314.

**BEFORE APPLYING ANY RESTRAINT:** Follow your hospital's restraint policies and procedures that are compliant with USA CMS guidelines and state laws or other governing agencies outside the USA.

**⚠️ ADDITIONAL WARNINGS:**

1. ALWAYS monitor patient per facility policy. Improper application or use of any restraint may result in serious injury or death. Be aware that constant monitoring may be required for:
  - Aggressive or agitated patients; and
  - Patients in a supine position; ensure that being in this position would not cause any complications for those who are at risk of vomiting while being restrained.
2. NEVER alter or repair this product. ALWAYS inspect before each use: Check for broken stitches or parts; torn, cut or frayed material; or locks, buckles, or hook-and-loop fasteners that do not hold securely. DO NOT use soiled or damaged products. Doing so may result in serious injury or death. Dispose of damaged products per facility policy for BIOHAZARDOUS material.
3. ALWAYS secure straps to a part of the frame that moves with the patient (not foot/head board or siderails), out of the patients reach, using quick-release ties (see drawing below) or buckles. These allow easy release in the event of an accident or fire.
 

Test to make sure straps cannot tighten, loosen, or slip and create excess slack. If this occurs, the patient may slide off the chair or bed, increasing the risk of serious injury or suffocation. Restraint release is an important part of facility fire and disaster drills. Straps can be cut with scissors in an emergency.
4. NEVER secure restraint strap to side rail or head/foot board.
5. NEVER use Posey products on toilets, or on any chair or furniture that does not allow proper application as directed in the Application Instructions. DO NOT use at home.
6. NEVER expose this product to open flame, fire, smoking materials, or high heat sources. Some products may melt or ignite and burn. The facility smoking/no smoking policy should be strictly enforced.
7. NEVER use a Posey product as a seat belt in a moving vehicle. Posey products are not designed to withstand the force of a crash or sudden stop.
8. Test Zippers or hook-and-loop fasteners before each use. DISCARD device if it does not fasten securely.

**LAUNDERING INSTRUCTIONS** (if applicable):

- Fasten all buckles and locks to reduce risk of damage during wash and dry cycles. DO NOT put buckles or locks through extractors. For maximum life, launder in a laundry bag.
- Stainless Steel Locks: Posey recommends that locks be treated with a dryfilm lubricant (silicone) after each cleaning. For best results, spray a small amount into the lock. Work the action of the lock with a key several times to ensure that all the surfaces are well lubricated.
- Before laundering, zip up and turn the product inside out to protect zipper.
- Hook-and-loop fasteners may collect lint after repeated use or laundering, reducing grip strength. Fasten the "hook" to the "loop" before laundering to help prevent lint buildup. As needed, use a stiff-bristle brush to remove lint from the "hook" side.

• These products, other than foam products, can be machine washed under CDC guidelines for material soiled with blood or bodily fluid.

- For non-contaminated material, use lower temperature wash and dry cycles to extend product life.

• For foam products:


- Make sure products are completely dry before use.

**STORAGE AND HANDLING:**

- This device is designed for use in normal indoor environments.
- This device may be stored in ambient warehouse temperatures at normal humidity levels. Avoid excess moisture or high humidity that may damage product materials.

**HOW TO TIE THE POSEY QUICK-RELEASE TIE**

1. Wrap strap once around the frame that moves with the patient (do not attach to side rail or head/footboard), leaving at least an 8" (20 cm) tail. Fold the loose end in half to create a loop and cross it over the other end.
2. Insert the folded strap where the straps cross over each other, as if tying a shoelace. Pull on the loop to tighten.
3. Fold the loose end in half to create a second loop.
4. Insert the second loop into the first loop.
5. Pull on the loop to tighten. Test to make sure strap is secure and will not slide in any direction.
6. Repeat on other side. Practice quick-release ties to ensure the knot releases with one pull on the loose end of the strap.

EN

Report any serious incident related to device use to TIDI Products and the Member State Competent Authority. | Translations of this document in French, Spanish, Dutch, German, Italian, and Portuguese are located at [www.tidiproducts.com/IFU](http://www.tidiproducts.com/IFU).

**DISPOSAL:** Dispose of product according to local or facility guidelines if product is worn or damaged.

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© 2025 TIDI Products, LLC. Alle rechten voorbehouden. Posey is een gedeponeerde handelsmerk van TIDI Products, LLC.
TIDI Products, LLC • 570 Enterprise Drive, Neenah, WI 54956 EUA
MDSS GmbH