INSTRUCTIONS

EN

Intended for use in the general population, for adhesive securing of the catheter hub, tube or line to the patient's skin.

- Select the correct Grip-Lok for your application. Follow your hospital's policies and procedures for dressing application, removal and assessment.

REF	For use with Winged Catheters	Includes Dressing Kit
3300MIV	✓	
3300MART		
3300MIV-KIT-2	✓	✓
3300MART-KIT-2		✓

2. Select the area for the placement of the Grip-Lok. Note: The catheter should be already inserted in the patient with Luer connection attached.

3. Prepare the skin according to the standard hospital protocol for dressing application. Hair removal may be required on some patients for better adhesion.

4. Open the Grip-Lok and fold back the top liner to expose center section of adhesive.

5. Slide the Grip-Lok under the catheter and the Luer connector centering in the exposed adhesive area.

6. Hold the Grip-Lok and catheter in position and remove the folded liner.

7. Hold the Grip-Lok in position and remove one side of the bottom liner then remove the other side of the bottom liner and gently press adhesive to skin.

8. Wrap top adhesive tab around the catheter pressing to conform to the shape of the catheter.

9. Wrap the remaining top adhesive tab around the Luer connector pressing to conform to the shape of the connector.

Direction: For additional skin protection, use skin prep pad prior to application of stabilization device.

Note: Use of an alcohol swab may help removal of securement device from the skin.

⚠️ If you experience severe redness, itching, swelling, or irritation of the skin, consult your physician as this may be a sign of an allergic reaction.

⚠️ Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

⚠️ Replace securement device if soiled or saturated in fluid or if device shows signs of wear or damage.

Note: If being used in a medical facility, dispose according to medical facility and municipal guidelines.

🚫 Do not use if the product's sterile barrier system or its packaging is compromised.

Report any serious incident related to device use to TIDI Products and the Member State competent authority

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HYPALLERGEN

TIDI Products, LLC
570 Enterprise Drive
Neenah, WI 54956, USA
Phone: +1.920.751.4300
tidiproducts.com/securement

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UK Responsible Person:
Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

Grip-Lok is a registered trademark of TIDI Products, LLC.
May be patented: see www.tidiproducts.com/virtual-patent-marking

2023-08-22

3300MIV

3300MIV-KIT-2

3300MART

3300MART-KIT-2

Grip-Lok®

Medium Peripheral IV and Arterial Catheter Securement Device

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