INSTRUCTIONS

EN

Intended for use in the general population, for adhesive securement of the catheter hub, tube or line to the patient's skin. English

1. Select the correct size Grip-Lok for your application.

REF	Size	Line/Tubing Diameter
3600PFC	Small	8 – 30 Fr
3400LFC	Large	8 – 30 Fr

2. Select the area for the placement of the Grip-Lok. Note: The catheter should already be inserted.
3. Prepare the skin according to the standard hospital protocol for dressing application. Skin prep or hair removal may be required on some patients for better adhesion.
4. Hold the Grip-Lok in position while removing the paper backing from one side, then the other, to secure on the skin surface.
5. Open and separate the top fabric straps.
6. Place the drainage lumen between the foam strip and the back of the fabric straps.
7. Position the catheter so that the Y-point is centered at the keyhole in the top fabric straps.
8. Secure one of the fabric straps over the drainage lumen and the other over the catheter below the Y-point.
9. Adjust the securement of the catheter by applying the strap attachments tighter or looser as needed.

Direction: For added skin protection, use skin prep pad prior to application of stabilization device.

Note: Use of an alcohol swab may help removal of securement device from the skin.

- ⚠ If you experience severe redness, itching, swelling, or irritation of the skin, consult your physician as this may be a sign of an allergic reaction.
- ⚠ Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.
- ⚠ Replace securement device if soiled or saturated in fluid or if device shows signs of wear or damage.
- ⚠ **Note:** If being used in a medical facility, dispose according to medical facility and municipal guidelines.
- ⊗ Do not use if the product's sterile barrier system or its packaging is compromised. Report any serious incident related to device use to TIDI Products and the Member State competent authority.

HYPOALLERGENIC

TIDI® Products, LLC
570 Enterprise Drive
Neenah, WI 54956 USA
Phone: +1 920 751 4300
tidiproducts.com/securement

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UK CA 0086

UK Responsible Person:
Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

Grip-Lok is a registered trademark of TIDI Products, LLC.

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Grip-Lok er et registrert varemærke tilhørende TIDI Products, LLC.