INSTRUCTIONS

EN

English

Intended for use in the general population.

1. **Compare** PICC/CVC wing shape to foam cavity in the Grip-Lok® to confirm compatibility.

REF	FOR
3304MCS-BA	Bard® PICC and CVC
3306MCS-NA	Navilyst Medical Xcela® PICC
3308MCS-MC	Medcomp® Dialysis Catheter
3303MCS-TA	Arrow® PICC/CVC

2. **Select** the area for the placement of the Grip-Lok. Note: The catheter should be already inserted in the patient.
3. **Prepare** the skin according to the standard hospital protocol for dressing application. Hair removal may be required on some patients for better adhesion.
4. **Open** the top flap, remove the interior liner and hold the flap open.
5. **Slide** the Grip-Lok under the catheter and adhere in the exposed adhesive cut-out space provided.
6. **Secure** the top flap over the catheter.
7. **Hold** the Grip-Lok in position while removing the paper backing from one side, then the other, to secure on the skin surface.

Direction: For additional skin protection, use skin prep pad prior to application of stabilization device.

Note: Use of an alcohol swab may help removal of securement device from the skin.

- ⚠ If you experience severe redness, itching, swelling, or irritation of the skin, consult your physician as this may be a sign of an allergic reaction.
- ⚠ Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.
- ⚠ Replace securement device if soiled or saturated in fluid or if device shows signs of wear or damage.
- Note:** If being used in a medical facility, dispose according to medical facility and municipal guidelines.
- ⊘ Do not use if the product's sterile barrier system or its packaging is compromised. Report any serious incident related to device use to TIDI Products and the Member State competent authority.

HYPOALLERGENIC

TIDI® Products, LLC
570 Enterprise Drive
Neenah, WI 54956 USA
Phone: +1 920 751 4300
tidiproducts.com/securement

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UK CA 0086

UK Responsible Person:
Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom


Grip-Lok is a registered trademark of TIDI Products, LLC.
May be patented: see www.tidiproducts.com/virtual-patent-marking

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