INSTRUCTIONS

EN

Intended for use in the general population.

- Select the area for the placement of the Grip-Lok. Note: The catheter should be already inserted in the patient.
- Prepare the skin according to the standard hospital protocol for dressing application. Skin prep or hair removal may be required on some patients for better adhesion.
- Open the top flap and hold the flap open.
- Position the Grip-Lok with the catheter in the groove at edge of the foam section.
- Hold the Grip-Lok in position while removing the paper backing from one side, then the other, to secure on the skin surface.
- Remove the interior liner and place the catheter over the foam and onto the center of the exposed adhesive surface.
- Secure the top flap over the catheter and apply gentle pressure to the center fabric area over the catheter adhesive.

Direction: for added skin protection, skin prep may be used prior to application of securement device.

Note: Use of an alcohol swab may help removal of securement device from the skin.

⚠ If you experience severe redness, itching, swelling, or irritation of the skin, consult your physician as this may be a sign of an allergic reaction.

⚠ Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

⚠ Replace securement device if soiled or saturated in fluid or if device shows signs of wear or damage.

Note: If being used in a medical facility, dispose according to medical facility and municipal guidelines.

⚠ Do not use if the product's sterile barrier system or its packaging is compromised.

Report any serious incident related to device use to TIDI Products and the Member State competent authority.

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TIDI
Support Caregivers.
Protect Patients.

TIDI Products, LLC
570 Enterprise Drive
Neenah, WI 54956, USA
Phone: +1.920.751.4300
tidiproducts.com/securement

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UK Responsible Person:
Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

STERILE R

Grip-Lok is a registered trademark of TIDI Products, LLC.
May be patented: see www.tidiproducts.com/virtual-patent-marking

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