INSTRUCTIONS

EN

English

Intended for use in the general population.

- Select the area for the placement of the Grip-Lok. Note: The feeding tube should be already inserted in the patient.
- Prepare the skin according to the standard hospital protocol for dressing application. Hair removal may be required on some patients for better adhesion.
- Remove the bottom liner by holding fabric section and pulling-off liner. Do not remove the top liner at this time.
- Place the Grip-Lok on the site by continuing to hold by the fabric and press onto the skin with gentle pressure. The top liner should now be removed.
- Position the feeding tube in the center of the top exposed adhesive strip between the fabric and the white hook section.
- Secure the feeding tube by folding over the fabric section until it meets with the white hook portion and apply gentle pressure to the back of the entire fabric area.

Note: to remove or adjust secured tube, hold hydrocolloid adhesive down to skin while opening top fabric section.

Note: the Grip-Lok 2100NGH may be cut in half prior to removal of liners to reduce the size.

Direction: For additional skin protection, use skin prep pad prior to application of stabilization device.

Note: Use of an alcohol swab may help removal of securement device from the skin.

- ⚠ If you experience severe redness, itching, swelling, or irritation of the skin, consult your physician as this may be a sign of an allergic reaction.
- ⚠ Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.
- ⚠ Replace securement device if soiled or saturated in fluid or if device shows signs of wear or damage.
- Note:** If being used in a medical facility, dispose according to medical facility and municipal guidelines.
- ⚠ Do not use if the product's sterile barrier system or its packaging is compromised. Report any serious incident related to device use to TIDI Products and the Member State competent authority.

TIDI® Products, LLC
570 Enterprise Drive
Neenah, WI 54956 USA
Phone: +1 920 751 4300
tidiproducts.com/securement

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UK Responsible Person:
Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom


Grip-Lok is a registered trademark of TIDI Products, LLC.
May be patented: see www.tidiproducts.com/virtual-patent-marking

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Grip-Lok je registrovaná ochranná známka společnosti TIDI Products, LLC.