

Make sure table top is clean and dry before applying drape.  
 Make sure there are no rips or tears.  
 Make sure holes in drape expose hook & loop Velcro to enable pads to be securely fastened.




**REF** BF508 x1

## STERIS® Surgical Sliding Tabletop Cover (Disposable)

**To Be Used With**  
 STERIS® 4000- & 5000-Series General Surgical Tables

**Intended Use:**  
 To cover and protect equipment.

⚠ Do not use if product is torn or punctured.



**To Re-Order, Please Contact:**

 STERIS Corporation  
 5960 Heisley Road, Mentor, OH 44060  
 440-354-2600 • 800-548-4873  
 www.steris.com

 TIDI Products, LLC  
 570 Enterprise Drive  
 Neenah, WI 54956, USA  
 1.800.521.1314  
 + 1.920.751.4300  
 www.tidiproducs.com

