

A New Generation Safety Scalpel

PenBlade is a self-retracting, safety scalpel designed for the physician. Utilizing an ergonomic grip and a high-quality scalpel blade, PenBlade makes precise surgical incisions with minimal effort. The safety mechanism is activated without changing the grip, creating a truly safe solution. PenBlade is safety without compromise.



“ PenBlade is the first safety scalpel suitable for use in the operating room. ”

- General Surgeon

“ It's perfectly sharp, easy to use, and much safer than a regular scalpel. ”

- Pediatric Surgeon

PenBlade Safety Advantages

Feature	Benefit
Self-Retracting Blade	For maximum safety the blade completely retracts with a push of a button
One-Handed Activation	The blade is easily deployed with one hand using a pen-like action
Front Release Button	The safety mechanism can be activated without looking, which allows clinicians to maintain focus on the patient
Integrated Trimming Groove	Clinicians can use the retracted blade for trimming sutures and PICCS
Effortless Safety	Minimizes the motion and effort required to operate the device safely
Rigid and Ergonomic	PenBlade is sturdy and easy to hold
Audible Click	An audible click indicates that the blade is fully protected and safe without looking



Part #	Blade Description	UOM	#/UOM
PB-M-10-CAS	#10 Sterile Safety Scalpels	Case	100
PB-M-11P-CAS	#11P Sterile Safety Scalpels	Case	100
PB-M-15-CAS	#15 Sterile Safety Scalpels	Case	100

GPO Contract Details

HealthTrust Purchasing Group: HPG-19318
 Intalere Healthcare: VH11573
 Partners Cooperative, Inc.: PA-001-212016
 Premier, Inc.: PP-OR-1665
 Vizient, Inc.: MS4610

+1 920.751.4300
www.tidiproductions.com

