


# Posey® Alarm Canister Bracket

## Application Instructions



### Application Instructions

1. Install 8208CB bracket in canister wall bracket.
2. Tighten set screw using a Phillips screwdriver until the 8208CB bracket is not removable (do not overtighten).
3. Install Posey alarm onto 8208CB bracket.




Set Screw



Alarm Canister Bracket

REF 8208CB Alarm Canister Bracket



TIDI Products, LLC • 570 Enterprise Drive, Neenah, WI 54956 USA  
Phone: 1.800.521.1314 • International: +1.920.751.4036 • [www.tidiproducts.com](http://www.tidiproducts.com)

© 2024 TIDI Products, LLC. All rights reserved. Posey is a registered trademark of TIDI Products, LLC.

M6164 REV B 03-20-2024


# Posey® Alarm Canister Bracket

## Application Instructions



### Application Instructions

1. Install 8208CB bracket in canister wall bracket.
2. Tighten set screw using a Phillips screwdriver until the 8208CB bracket is not removable (do not overtighten).
3. Install Posey alarm onto 8208CB bracket.




Set Screw



Alarm Canister Bracket

REF 8208CB Alarm Canister Bracket



TIDI Products, LLC • 570 Enterprise Drive, Neenah, WI 54956 USA  
Phone: 1.800.521.1314 • International: +1.920.751.4036 • [www.tidiproducts.com](http://www.tidiproducts.com)

© 2024 TIDI Products, LLC. All rights reserved. Posey is a registered trademark of TIDI Products, LLC.

M6164 REV B 03-20-2024