

Make sure table top is clean and dry before applying drape.
 Make sure there are no rips or tears.
 Make sure holes in drape expose mushroom cap fasteners to enable pads to be securely fastened.




REF BF549 x1

STERIS® Surgical Sliding Table Drape For Velcro®-less Tabletops (Disposable)

To Be Used With
 STERIS® 4000- & 5000-Series General Surgical Tables


Intended Use:
 To cover and protect equipment.

 Do not use if product is torn or punctured.



To Re-Order. Please Contact:

 STERIS Corporation
 5960 Heisley Road, Mentor, OH 44060
 440-354-2600 • 800-548-4873
 www.steris.com

 TIDI Products, LLC
 570 Enterprise Drive
 Neenah, WI 54956, USA
 1.800.521.1314
 + 1.920.751.4300
 www.tidiprotucts.com

