

WIRELESS

FALL ALARMS AND SENSORS

Support Caregivers. Protect Patients.

Why Go Wireless?



Patient Safety

Fewer cords in the room, fewer tripping hazards for the patient



Staff Efficiency

Ease of setup, ease of use, improved workflow, flexible alarm mounting



Financial Impact*

The cost of adding a wireless fall alarm or sensor is less than the cost of a single patient fall

Commonly Asked Questions

Do Posey alarms and sensors collect data?

No, Posey alarms and sensors do not collect any data.

Do Posey alarms and sensors use a Wi-Fi® wireless network?

No, Posey alarms and sensors do not use the hospital's wireless network, e.g., Wi-Fi® brand or otherwise, they connect to one another via Bluetooth® wireless communication protocol.



Trusted Technology. Proven Protection.

Wireless Toilet Sensor

Helps prevent unassisted toilet exits and boosts fall prevention protocols.

Wireless Chair Sensor Pad

Allows caregivers to move the chair anywhere in the room, improving both room layout and nurse workflow.

Wireless On Cue® PRO Alarm

Can be placed anywhere in the room; easily pairs with the Posey Wireless Chair Sensor Pad, Wireless Toilet Sensor, and Wireless Nurse Call Adapter.

Wireless Nurse Call Adapter

Pairs with the Posey On Cue PRO Alarm quickly and simply so that when the alarm sounds, the nurse call system will be wirelessly notified.

* Note: Direct costs of noninjurious and injurious inpatient falls average as much as \$35,365 and \$36,776 per fall, respectively. [Dykes PC, Curtin-Bowen M, Lipsitz S, et al. Cost of inpatient falls and cost-benefit analysis of implementation of an evidence-based fall prevention program. JAMA Health Forum. 2023;4(1):e225125. doi:10.1001/jamahealthforum.2022.5125]

© 2026 TIDI Products, LLC. All rights reserved.

Wi-Fi® is a registered trademark of and owned by the Wi-Fi Alliance®, relates to the IEEE 802.11 family of wireless networking standards, and is not associated or affiliated with TIDI Products, LLC.

Bluetooth® is a registered trademark of and owned by Bluetooth SIG, Inc., and while used in TIDI® products is not associated or affiliated with TIDI Products, LLC.

On Cue®, Posey®, and TIDI® are registered trademarks of and owned by TIDI Products, LLC.

Ready to Learn More?

Contact a Sales Rep

